

Title of the invention : A uric acid diagnostic kit and a method of detecting uric acid

International classification

:A61P0019060000,
G01N0033543000,
A61K0009000000,
C12Q0001620000,
C12Q0001688600
:NA
:NA
:NA
:NA
:NA
:NA
:NA
:NA
:NA

Priority Document No

Priority Date

Name of priority country

International Application No

Filing Date

International Publication No

Patent of Addition to Application Number

Filing Date

Divisional to Application Number

Filing Date

(71)Name of Applicant :

1)Misar Pranaya Laxman

Address of Applicant :N.N. Sattha college of Pharmacy, Ahmednagar-414001, Maharashtra, India Maharashtra India

2)Rachana Hajare

3)Payal Navsupe

4)Kuskar Sharmilla

5)Shrutika Dhananjay Dalvi

(72)Name of Inventor :

1)Shreya Dhadge


2)Anushka Pandit

3)Anant Bital

Abstract :

TRACT A uric acid diagnostic kit and a method of detecting uric acid The present invention is a uric acid diagnostic kit and a method of making thereof. The kit ensures the uric acid level effectively and establishes a robust and scalable method for immobilizing uricase enzyme on whatsmann filter paper, ensuring stability and effectiveness in salivary uric acid detection. The present invention pertains to a ground breaking uric acid determination kit having a unique electrochemical uric acid sensor which ensures swift and precise results. The kit is a comprehensive, efficient, and economical device, poised to revolutionize at-home salivary uric acid monitoring and the method is not only champions non-invasiveness, making it inherently patient-friendly, affordable and reliable solution for uric acid detection. Figure


Preparation of membrane with 0.1% Polystyrene



Addition of Uricase + GSA on membrane




Preparation of membrane with 0.1% Polystyrene + GSA



Immobilization of Uricase on membrane (Kimmco) (Kimmco) (Kimmco) as per protocol described in literature

Notes: The protocol is improved with reference to additional literature research



Packaging of uricase immobilized membrane



Membrane washing and incubation with 0.1M Glycine

FIGURE 1